English

Quinn's Advantage™ Sperm Freeze

For laboratory procedures only; other uses must be qualified by the end user.

Product Description	REF Number	Unit Size
Quinn's Advantage™ Sperm Freeze	ART-8022	6 x 12 mL

INTENDED USE

This product is intended for freezing sperm.

PRODUCT DESCRIPTION

Quinn's Advantage™ Sperm Freeze is a HEPES-buffered salt solution containing 10 mg/mL Human Serum Albumin, glycerol as the cryoprotective agent (CPA), and gentamicin as an antibiotic. Glycerol has been used as a cryoprotectant for human spermatozoa for nearly 50 years (Larson et al, 1997; Quinn, 1993). Similar glycerol/HEPES-buffered formulations demonstrated higher postthaw motility, longevity, percent survival, and ability to penetrate cervical mucus than conventional methods (Larson et al, 1997).

This product contains 10mg/mL of genta-micin, an aminoglycoside antibiotic.

PRECAUTIONS AND WARNINGS

Do not use medium that shows evidence of particulate matter.

This product contains 10 mg/mL human serum albumin.

Caution: All blood products should be treated as potentially infectious. Source material from which this product was derived was found negative when testing for antibodies to HIV-1/ HIV-2, HCV and non-reactive for HBsAg, HCV RNA and HIV-1 RNA. No known test methods can offer assurances that products derived from human blood will not transmit infectious agents. Donors of the source material have been screened for Creutzfeldt-Jacobs disease (CJD). Based on effective donor screening and product manufacturing processes, it carries an extremely remote risk for transmission of viral diseases. A theoretical risk for transmission of CJD is also considered extremely remote. No cases of transmission of viral diseases or CJD have ever been identified for albumin.

Standard measures to prevent infections resulting from the use of medicinal products prepared from human blood or plasma include selection of donors, screening of individual donors and plasma pools for specific markers of infection and the inclusion of effective manufacturing steps for the inactivation/removal of viruses. Despite this, when medicinal products prepared from human blood or plasma are administered, the possibility of transmitting infective agents cannot be totally excluded. This also applies to unknown or emerging viruses and other pathogens. There are no reports of proven virus transmissions with albumin manufactured to European Pharmacopoeia specifications by established processes.

Single use: To avoid problems with contamination, handle using aseptic techniques and discard any excess product that remains in the bottle or vial after procedure is completed.

Reproductive media products are intended for single use only. Re-use of reproductive media may result in using a product past its labeled

expiration date or increase the risk of microbial contamination in a subsequent procedure if the practitioner fails to utilize adequate aseptic techniques.

Use of expired or microbial contaminated product may result in suboptimal conditions to promote fertilization and/or embryo quality during in-vitro culture. These conditions may result in the failure of the embryo to develop properly or to implant, potentially leading to a failed assisted reproductive procedure.

Caution: U.S. Federal law restricts this device to sale by or on the order of a physician (or properly licensed practitioner).

This product contains the antibiotic genta-micin sulfate. Appropriate precautions should be taken to ensure that the patient is not sensitized to this antibiotic.

QUALITY ASSURANCE

One-cell MEA tested and passed with 80% or greater blastocyst. USP Endotoxin tested and passed with <1 EU/mL.

A Certificate of Analysis is available for this product.

DIRECTIONS FOR USE

1. Semen samples are allowed to liquefy in an incubator at 37 °C for 30 minutes after collection. Washed spermatozoa can be prepared with a density-gradient separation medium (PureCeption™ 40% [REF #2040] or 80% [REF #2080]) and the final washed sperm preparation resuspended in Quinn's™ Sperm Washing Medium (REF #1005/1006). Alternatively, the liquefied semen sample can be cryopreserved directly.
2. Bring all reagents to room temperature (20 °C to 25 °C) before use.

FREEZING

1. Add one volume of Quinn's Advantage™ Sperm Freeze dropwise over a 30 second period to one volume of liquefied semen or washed spermatozoa solution, with continual mixing after each drop of Quinn's Advantage™ Sperm Freeze is added. It is important to add the medium dropwise over the 30 second period and to mix the Quinn's Advantage™ Sperm Freeze thoroughly after each drop is added to allow for adequate equilibration of the freezing medium with the sperm cells. Allow the Quinn's Advantage™ Sperm Freeze/ sperm mixture to equilibrate for 3 minutes after all the Quinn's Advantage™ Sperm Freeze has been added.

2. Place the Quinn's Advantage™ Sperm Freeze/sperm mixture into straws or vials and cool at 0.5 °C/minute from 25 °C to -5 °C. Hold the cryocontainer at -5 °C for 3 minutes. Seed the cryocontainers manually by touching them for about one second with forceps precooled in liquid nitrogen. Hold the cryocontainers at -5 °C for another 7 minutes. Cool the cryocontainers at 10 °C/ minute from -5 °C to -80 °C. Plunge the cryo-ntainers in liquid nitrogen and then transfer them to storage canes.

3. Alternatively, suspend the vials/straws on aluminum canes and immerse in a container filled with approximately 600 mL of water at room temperature. Place the container of water holding the samples in a refrigerator at 4 °C for 30 to 90 minutes.

4. Transfer quickly to liquid nitrogen vapor at the top of a liquid nitrogen storage tank and leave for 30 to 45 minutes. Vials should be suspended about 10 to 20 cm above the surface of the liquid nitrogen. Straws should be

laid horizontally at a similar height. Quickly transfer the vials/straws to final storage on labeled canes in liquid nitrogen.

5. The next day, or several hours later, thaw a test vial/straw and record all results on an appropriate report form.

THAWING, DILUTION, AND WASHING

1. Thaw straws by placing them on the bench-top (22 °C). Cryovials need to be agitated in a water bath at 30 °C to 35 °C.

2. Transfer the thawed sperm suspension to a culture tube of adequate volume and then slowly add dropwise 10 volumes of Quinn's Sperm Washing Medium (REF #1005/1006) to the thawed sperm suspension over a 30 second period with adequate mixing to ensure complete dilution of the Quinn's Advantage™ Sperm Freeze.

Motile spermatozoa are then recovered from the thawed, diluted suspension by density-gradient centrifugation and washing.

Each laboratory should make its own determination of which medium to use for each particular procedure.

Information on specific aspects of IVF, embryo culture, and cryopreservation is available in our Product Catalog.

STORAGE INSTRUCTIONS AND STABILITY

Unopened containers must be stored at 2 °C to 8 °C upon receipt. Warm to incubator (37 °C) temperature prior to use. Quinn's Advantage™ Sperm Freeze should be tightly capped when used in a CO₂ incubator to avoid downward drifts of pH. Do not expose to temperatures greater than 39 °C. The product is stable until the expiration date shown on the label.

A. Remove desired volume of product using aseptic procedures.

B. Once removed, do not return any volume of product to the original container.

C. Do not use if the product becomes discolored, cloudy, turbid, or shows any evidence of microbial contamination.

RELATED PRODUCTS

ART-2040 PureCeption 40% Upper Phase

ART-2080 PureCeption 80% Lower Phase

ART-1005/1006 Quinn's™ Sperm Washing Medium

SAGE In Vitro Fertilization™ has a full line of products for the Reproductive Medicine Specialist. Please call or write for specific information or to receive a copy of our current catalog. For technical questions, or to reach our Customer Service Department, call the SAGE™ Support Line.

REFERENCES

Larson JM, McKinney KA, Mixon BA, et al. An intraterrine insemination-ready cryopreservation method compared with sperm recovery after conventional freezing and post-thaw processing. Fertil Steril. 1997;68:143-148.

Quinn P. Cryopreservation. In: Assisted Reproductive Technologies. Boston, Mass: Blackwell Scientific Publications; 1993:89-107.

Call the SAGE™ SUPPORT LINE: In the U.S.: (800) 243-2974 International: (203) 601-9818

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Références : Larson JM, McKinney KA, Mixon BA, et al. An intraterrine insemination-ready cryopreservation method compared with sperm recovery after conventional freezing and post-thaw processing. Fertil Steril. 1997;68:143-148. Quinn P. Cryopreservation. In: Assisted Reproductive Technologies. Boston, Mass: Blackwell Scientific Publications; 1993:89-107.

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EXPLANATION OF SYMBOLS

	Catalog Number
	Batch Number
	Use By (year, month, day)
	Do Not Reuse
	Temperature Limitation
	Aseptic Technique Sterilization Membrane Filtered (SAL 10 ⁻³)
	ATTENTION: See instructions for use.
	Authorized Representative in the European Community.
	Product conforms to the Medical Device Directive 93/42/EEC
	Manufacturer

RX ONLY U.S. Federal law restricts this device to sale by or on the order of a physician (or properly licensed practitioner).

SAGE™ In Vitro Fertilization, Inc. a CooperSurgical Company
95 Corporate Drive
Trumbull, CT 06611 USA

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Bibliografia Larson JM, McKinney KA, Mixon BA, et al. An intraterrine insemination-ready cryopreservation method compared with sperm recovery after conventional freezing and post-thaw processing. Fertil Steril. 1997;68:143-148. Quinn P. Cryo-preservation. In: Assisted Reproductive Technologies. Boston, Mass: Blackwell Scientific Publications; 1993:89-107.

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In Vitro Fertilization
a CooperSurgical Company

SAGE™ In Vitro Fertilization, Inc.
a CooperSurgical Company
95 Corporate Drive
Trumbull, CT 06611 USA

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REFERENCE

Larson JM, McKinney KA, Mixon BA, et al. An intrauterine insemination-ready cryopreservation method compared with sperm recovery after conventional freezing and post-thaw processing. Fertil Steril. 1997;68:143-148.

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	In Vitro Fertilization в CooperSurgical Company
	SAGE™ In Vitro Fertilization, Inc. а CooperSurgical Company 95 Corporate Drive Trumbull, CT 06611 USA
	ORIGIO a/s Knardrupvej 2 2760 Mдløv Denmark www.fertility.coopersurgical.com
Customer Service:	E-mail: sales@coopersurgical.com Tel: +45 46 79 02 00

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