## English

### Quinn's Advantage™ Medium with HEPES

For laboratory procedures only; other uses must be qualified by the end user.

Product Description	REF Number	Unit Size
Quinn's Advantage™ Medium w/HEPES	ART-1023	100 mL
Quinn's Advantage™ Medium w/HEPES	ART-1024	500 mL

#### INTENDED USE

This product is intended for in vitro procedures involving manipulation of gametes and embryos not requiring the use of a CO<sub>2</sub> incubator.

#### PRODUCT DESCRIPTION

This medium is a modification of Human Tubal Fluid (HTF) that was first described by Quinn and colleagues (Fertil Steril. 1984;41:202, 1985;44:493.). Modifications to the original formulation include low phosphate and the addition of citrate, EDTA, selected nonessential amino acids and taurine, alanyl-glutamine as a stable source of glutamine, increased levels of magnesium, and the presence of lactate in the form of calcium lactate. This product contains 10 mg/L of gentamicin, an aminoglycoside antibiotic.

#### PRECAUTIONS AND WARNINGS

Do not use medium that shows evidence of particulate matter, cloudiness, or is not rose colored.

**Single Use:** To avoid problems with contamination, handle using aseptic techniques and discard any excess product that remains in the bottle or vial after procedure is completed.

Reproductive media products are intended for single use only. Re-use of reproductive media may result in using a product past its labeled expiration date or increase the risk of microbial contamination in a subsequent procedure if the practitioner fails to utilize adequate aseptic techniques.

Use of expired or microbial contaminated product may result in suboptimal conditions to promote fertilization and/or embryo quality during in-vitro culture. These conditions may result in the failure of the embryo to develop properly or to implant, potentially leading to a failed assisted reproductive procedure.

**Note:** Embryo is considered a general term. More precisely, SAGE™ considers the period of time initiating when a single diploid cell results from the fusion of male and female genome resulting in zygote formation with subsequent development from repeated mitotic divisions forming a solid mass or morula (typically day 4-5) and after which a fluid-filled cavity develops resulting in blastocyst formation (typically day 5-6) ending with embryo implantation that begins the end of the first week and is completed by the end of the second week post conception.

**Caution:** U.S. Federal law restricts this device to sale by or on the order of a physician (or properly licensed practitioner).

This product contains the antibiotic gentamicin sulfate. Appropriate precautions should be taken to ensure that the patient is not sensitized to this antibiotic.

#### QUALITY ASSURANCE

One-cell MEA tested and passed with 80% or greater blastocyst. USP Endotoxin tested and passed with <1 EU/mL.

A Certificate of Analysis is available for this product.

#### DIRECTIONS FOR USE

This is the preferred medium for procedures done in an atmosphere of air not involving the use of a CO<sub>2</sub> incubator. Such procedures include oocyte recovery, gamete washing, micromanipulation, embryo transfer, and cryopreservation. The medium should be supplemented with an exogenous protein source such as Human Serum Albumin (REF #3001/3003) or Serum Protein Substitute (REF #3010/3011).

Each laboratory should make its own determination of which medium to use for each particular procedure.

Information on specific aspects of IVF, embryo culture, and cryopreservation is available in our Product Catalog.

#### STORAGE INSTRUCTIONS AND STABILITY

Store unopened containers refrigerated at 2 °C to 8 °C. Warm to incubator (37 °C) temperature prior to use. Quinn's Advantage™ Medium with HEPES should be tightly capped when used in a CO<sub>2</sub> incubator to avoid downward drifts of pH.

Do not freeze or expose to temperatures greater than 39 °C. The product is stable until the expiration date shown on the label.

- Remove desired volume of product using aseptic procedures.
- Once removed, do not return any volume of product to the original container.
- Do not use if the product becomes discolored, cloudy, turbid, or shows any evidence of microbial contamination.

#### RELATED PRODUCTS

ART-3001/3003 Human Serum Albumin 100 mg/mL in normal saline  
ART-3010/3011 Quinn's Advantage SPS Serum Protein Substitute  
SAGE In Vitro Fertilization™ has a full line of products for the Reproductive Medicine Specialist. Please call or write for specific information or to receive a copy of our current catalog. For technical questions, or to reach our Customer Service Department, call the SAGE™ Support Line.

Quinn's Advantage™ is a trademark of CooperSurgical, Inc.

**Call the SAGE™ SUPPORT LINE:  
In the U.S.: (800) 243-2974 International:  
(203) 601-9818**

#### EXPLANATION OF SYMBOLS

- Catalog Number
- Batch Number
- Use By (year, month, day)
- Do Not Reuse
- Temperature Limitation
- Aseptic Technique Sterilization Membrane Filtered (SAL 10<sup>-3</sup>)
- ATTENTION:** See instructions for use. Authorized Representative in the European Community.
- Product conforms to the Medical Device Directive 93/42/EEC
- Manufacturer
- U.S. Federal law restricts this device to sale by or on the order of a physician (or properly licensed practitioner).

In Vitro Fertilization  
a CooperSurgical Company  
SAGE™ In Vitro Fertilization, Inc.  
a CooperSurgical Company  
95 Corporate Drive  
Trumbull, CT 06611 USA

ORIGIO a/s  
Knardrupvej 2  
2760 Måløv  
Denmark  
www.fertility.coopersurgical.com

**Customer Service:**  
E-mail: sales@coopersurgical.com  
Tel: +45 46 79 02 00

5627-06-ver. 6: 2019.Oct.17

## Français

### Quinn's Advantage™ Medium with HEPES

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Quinn's Advantage™ é uma marca registada da CooperSurgical, Inc.