English

**SAGE™ Vitrification  
Warming Kit**

For laboratory procedures only; other uses must be qualified by the end user.

| Product Description            | REF Number |
|--------------------------------|------------|
| 1.0 M Sucrose Warming Solution | ART-8030-A |
| 0.5 M Sucrose Warming Solution | ART-8030-B |
| MOPS Solution                  | ART-8030-C |

**INTENDED USE**

The products are intended for the recovery of human embryos (pronuclear zygotes through day 3 cleavage stage embryos and blastocysts stage embryos) that have been vitrified using the SAGE™ IVF Vitrification Kit (Ref # ART-8025). This kit is designed to be used in conjunction with the SAGE™ Vitrification Kit (ART-8025) for ultra-rapid freezing of specimens.

**PRODUCT DESCRIPTION**

1.0 M Sucrose Warming Solution (ART-8030-A) is a MOPS buffered solution of modified HTF containing non-essential and essential amino acids, gentamicin sulfate (10 mg/L), 1.0 M sucrose and 12 mg/mL human serum albumin.

0.5 M Sucrose Warming Solution (ART-8030-B) is a MOPS buffered solution of modified HTF containing non-essential and essential amino acids, gentamicin sulfate (10 mg/L), 0.5 M sucrose and 12 mg/mL human serum albumin.

MOPS Solution (ART-8030-C) is a MOPS buffered solution of modified HTF containing

non-essential and essential amino acids, gentamicin sulfate (0.01 g/L) and 12 mg/mL human albumin.

This product contains 10 mg/L of gentamicin, an aminoglycoside antibiotic.

**MATERIALS PROVIDED IN THE VITRIFICATION WARMING KIT**

1 x 4 mL vial of 1.0 M Sucrose Warming Solution (REF # ART-8030-A)

1 x 4 mL vial of 0.5 M Sucrose Warming Solution (REF # ART-8030-B)

PRECAUTIONS AND WARNINGS

**Warning:** The long term safety of embryo(s) vitrification on children from this procedure is unknown.

**Caution:** The safety and effectiveness of vitrification has not been fully evaluated in human embryos that have not yet reached the cleavage stage. Data reported from cleavage-stage embryos vitrified on day 3 of development (Desai et al., 2007) and clinical pregnancies have been established from vitrified zygotes (Selman & El-Danoush, 2002).

Standard measures to prevent infections resulting from the use of medical products prepared from human blood or plasma include selection of donors, screening of individual donations and plasma pools for specific markers of infection and the inclusion of effective manufacturing steps for the inactivation/inactivation of viruses.

Despite these measures, some medical products prepared from human blood or plasma are administered, the possibility of transmitting infective agents cannot be totally excluded. This also applies to unknown or emerging viruses and other pathogens. There are no reports of proven virus transmissions with albumin manufactured to European Pharmacopoeia specifications by established processes.

A Certificate of Analysis is available for this product.

**Caution:** The user should read and understand the Directions for Use. Precautions and Warnings, and be trained in the correct procedure before using the Vitrification and Warming Kits for vitrification of human embryos.

The 1.0 M Sucrose Warming Solution (ART-8030-A) component in this kit contains 12 mg/mL human serum albumin.

The 0.5 M Sucrose Warming Solution (ART-8030-B) component in this kit contains 12 mg/mL human serum albumin.

The MOPS Solution Solution (ART-8030-C) component in this kit contains 12 mg/mL human serum albumin.

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