Embryo Replacement Catheter

en

INSTRUCTIONS FOR USE

fr

MODE D'EMPLOI

de

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it

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es

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pt

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nl

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sv

BRUKSANVISNING

da

BRUGSANVISNING

no

BRUKSANVISNING

fi

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el

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CE
0473

 Manufacturer:
Smiths Medical International Ltd.
Boundary Road, Hythe
Kent, CT21 6JL, UK
Tel: +44 (0)1233 722100

www.smiths-medical.com

CCN. 004/006/007A 07/14

Wallace® Embryo Replacement Catheters

INSTRUCTIONS FOR USE

Wallace® embryo replacement catheters are single-use devices provided sterile for the introduction of embryos into the uterine cavity following in vitro fertilisation.

These instructions contain important information for safe use of the product. Read the entire contents of these Instructions For Use, including Warnings and Cautions, before using this product. Failure to properly follow warnings, cautions and instructions could result in death or serious injury to the patient and/or clinician.

NOTE: DISTRIBUTE THIS INSERT SHEET TO ALL PRODUCT LOCATIONS.

WARNING:

- **THESE PROCEDURES SHOULD ONLY BE PERFORMED BY PERSONS HAVING ADEQUATE TRAINING AND FAMILIARITY WITH THESE TECHNIQUES. CONSULT MEDICAL LITERATURE REGARDING TECHNIQUES, COMPLICATIONS AND HAZARDS PRIOR TO PERFORMANCE OF THESE PROCEDURES.**
- **TO BE USED BY, OR UNDER THE DIRECTION OF, QUALIFIED PERSONS IN LINE WITH LOCAL GUIDELINES GOVERNING IN VITRO FERTILISATION, IF APPLICABLE.**
- **DISCARD IF PRODUCT OR PACKAGING IS DAMAGED.**

Description

These instructions are intended for use with sterile Wallace® embryo replacement catheters, product codes

Catheter type	Length	Stiffness	Code
Classic	18cm	Soft	1816
Classic	23cm	Soft	1816N
Classic	18cm	Intermediate	1816S
Classic	23cm	Intermediate	2316
Sure View**	18cm	Soft	CE118
Sure View**	23cm	Soft	CE123

*Sure View® – for enhanced visibility under ultrasound.

The devices consist of a flexible end opening inner catheter and a detachable outer sheath that is attached to the inner catheter by a Luer compatible adapter. The inner catheter is 16g with a lumen that is uniform throughout its length, with a series of 1cm graduations at the proximal end. The inner catheter protrudes from the outer sheath by 5cm and the outer sheath has a series of 1cm graduations at the distal end. The material of the Sure View® inner catheter is modified to contain small bubbles that improve ultrasound image quality.

Indications for use

Introduction of embryos into the uterine cavity.

Contraindications

The catheter should not be used:

- In the presence of chronic cervical infection
- In the presence of or after recent pelvic inflammatory disease
- For intra-fallopian procedures

Cautions

- Low resolution ultrasound equipment, or use of 'Classic' embryo catheters may compromise the visibility of the catheter.
- When performing transabdominal ultrasound, ensure that the patient has a full bladder.
- The patient should be fully prepared for embryo transfer before preparation and loading of the catheter commences.

Instructions for Use

Place patient in the lithotomy, dorsal, or left lateral position. The cervix should be exposed with a speculum and gently swabbed with cotton wool moistened with normal saline or medium.

When using transabdominal ultrasound follow normal procedures such that an optimal image of the uterine cavity is achieved.

Under sterile conditions, fill a 1ml or Hamilton syringe with medium and eject any air bubbles. Remove the tip protector and firmly attach the syringe to the Luer compatible adapter of the inner catheter. Flush the catheter through with medium, remove syringe, and refill with medium. Reconnect the syringe and expel medium until the 0.1ml mark is reached, with the catheter tip pointing upwards to facilitate the removal of air bubbles.

Place the end of the catheter carefully into the medium droplet containing the embryos and inject a small amount of medium to break the boundary of surface tension that may appear.

Aspirate the embryos into the catheter and hand to person performing embryo transfer.

With the hubs locked in position, advance the catheter so that the inner catheter passes through the external and internal os, into the mid-uterine cavity. It may be necessary to twist the catheter as it negotiates the cervical canal. When using ultrasound manipulate the probe to obtain the best image of the inner catheter.

WARNING

The catheter should never be forced against digitally felt resistance while inside the uterine cavity, as forcing the catheter may result in damage to the endometrial tissue and bleeding.

The outer sheath should not be advanced further than the internal os, and should certainly never enter the uterine cavity, as this may result in damage to the endometrial tissue and bleeding.

The inner catheter should not be pulled back further than the last graduation. Removing it from the outer sheath completely may result in damage to the endometrial tissue and bleeding.

If slight resistance is experienced, withdraw the catheter and advance the outer sheath until only the tip of the inner catheter is visible. Aligning the most distal marking on the base of the inner catheter with the hub of the outer sheath ensures a smooth radius protrudes. Using thumb and forefinger curve the outer sheath to complement the patient's anatomy and pass the catheter in this position to the internal os. In the rare event where the cervical canal cannot be negotiated by this method, stylets are available to assist insertion. Use code no 1816ST with 18cm catheters and 1816NST with 23cm catheters. Advance the inner catheter into the mid-uterine cavity.

Place tip of catheter approximately 1cm from the fundus. Twist the catheter to remove any mucus from the tip.

Expel the embryos.

Leave time for the uterine environment to equilibrate and slowly remove the catheter.

Hand catheter back to embryologist to check that no embryos remain. Remove the speculum and allow the patient to rest.

Dispose of the catheter in accordance with local medical hazardous waste practices.

The above are only suggested instructions for use. The final decision on the technique used is the responsibility of the clinician in charge.

*The products described are covered by one or more of the following: U.S. Patent No. 8,092,390 and counterpart foreign patent(s). Pat pending. Patent applied for.

Sure View, Wallace, Wallace design mark and Smiths Medical design mark are trademarks of Smiths Medical. The symbol ® indicates the trademark is registered in the U.S. Patent and Trademark Office and certain other countries. All other names and marks mentioned are the trade names, trademarks or service marks of their respective owners. © 2014 Smiths Medical. All rights reserved.

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Sure View, Wallace, designmerket Wallace og designmerket Smiths Medical er varemerker som tilhører Smiths Medical. Symbolet® angir at varemerket er registrert ved U.S. Patent and Trademark Office og i visse andre land. Alle andre navn og merker som er nevnt, er varemerker, handelsnavn og servicemerker som tilhører de respektive eiere.