Oocyte Recovery Sets

en

INSTRUCTIONS FOR USE

fr

MODE D'EMPLOI

de

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it

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es

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pt

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nl

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sv

BRUKSANVISNING

da

BRUGSANVISNING

no

BRUKSANVISNING

fi

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el

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CE
0473

 Manufacturer:
Smiths Medical International Ltd.
Boundary Road, Hythe
Kent, CT21 6JL, UK
Tel: +44 (0)1233 722100

www.smiths-medical.com

smiths medical

Wallace®

Oocyte Recovery Sets

INSTRUCTIONS FOR USE

These instructions apply to the use of the following Wallace® products:

Code	Description
ONS1633; ONS1633-500~750	16g 33cm needle with transfer tubing
ONS1733; ONS1733-500~750	17g 33cm needle with transfer tubing
ONS1833; ONS1833-500~750	18g 33cm needle with transfer tubing
ONS1633LL; ONS1633LL-500~750	16g 33cm needle, transfer tubing with flushing connector
ONS1733LL; ONS1733LL-500~750	17g 33cm needle, transfer tubing with flushing connector
ONS1833LL; ONS1833LL-500~750	18g 33cm needle, transfer tubing with flushing connector

These instructions contain important information for safe use of the product. Read the entire contents of these Instructions For Use, including Warnings and Cautions, before using this product. Failure to properly follow warnings, cautions and instructions could result in death or serious injury to the patient and/or clinician.

NOTE: DISTRIBUTE THIS INSERT SHEET TO ALL PRODUCT LOCATIONS.

WARNING:

- **THESE PROCEDURES SHOULD ONLY BE PERFORMED BY PERSONS HAVING ADEQUATE TRAINING AND FAMILIARITY WITH THESE TECHNIQUES. CONSULT MEDICAL LITERATURE REGARDING TECHNIQUES, COMPLICATIONS AND HAZARDS PRIOR TO PERFORMANCE OF THESE PROCEDURES.**
- **TO BE USED BY, OR UNDER THE DIRECTION OF, QUALIFIED PERSONS IN LINE WITH LOCAL GUIDELINES GOVERNING IN VITRO FERTILISATION, IF APPLICABLE.**
- **DISCARD IF PRODUCT OR PACKAGING IS DAMAGED.**

DESCRIPTION

Wallace® Oocyte Recovery Sets are single-use devices provided sterile for ultrasound-guided, transvaginal collection of oocytes from the ovarian follicles.

The devices consist of a single-lumen stainless steel needle with 2cm of echo marking at the distal end. Attached to the proximal end are a plastic hub and a tubing set. The tubing set consists of a transfer tube leading through a silicone bung and a vacuum connection tube also passing through the bung and terminating in a Luer compatible connector. Product codes with a LL suffix have an additional Luer connector on the transfer tube to enable flushing of the device.

INDICATIONS FOR USE

Ultrasound-guided, transvaginal collection of oocytes.

CONTRAINDICATIONS

The oocyte recovery set should not be used:

- In the presence of or after recent pelvic inflammatory disease.
- For intra-fallopian procedures.

CAUTIONS

Low resolution ultrasound equipment may compromise the visibility of the needle.

INSTRUCTIONS FOR USE

- 1 Oocyte recovery should be performed under either general anaesthetic or analgesia.
- 2 Place patient in the lithotomy position, prepare the vulva and vagina for full aseptic technique and cover the area in sterile drapes.
- 3 Cover the ultrasound transducer with a sterile sleeve with a small amount of coupling gel covering the end of the probe within the sleeve. Introduce into the vagina and scan the pelvis to establish the position of the uterus, the development of the endometrium and the number and position of follicles to be aspirated.

- 4 Remove the oocyte recovery set from its packaging and insert the needle into the needle guide of the ultrasound transducer, ensuring that there is smooth movement along its length.
- 5 Attach a syringe or vacuum pump to the Luer compatible connector on the vacuum line and a sterile collection tube to the silicone bung. The transfer tube can be pulled through the bung so that it reaches the bottom of the collection tube, preventing damage to the zona pellucida of the aspirated oocytes.

WARNING:

Test the system to ensure all connections are sound and flow can be achieved. Failure to test the system may result in failure to obtain the required sample and/or tissue damage.

- 6 Insert the ultrasound transducer and needle into the vagina and direct the probe towards the posterolateral aspect of the vaginal fornix. Manipulate the probe until the follicles appear in close proximity to the needle.

WARNING:

Failure to use ultrasound to locate the follicles may result in failure to obtain the required sample and/or tissue damage.

- 7 Advance the needle into the nearest follicle. If the ovaries are mobile apply pressure suprapubically and enter the follicles with a stabbing movement. When the tip of the needle can be seen in the follicle, aspirate the contents until empty using either the syringe or vacuum pump attached to the set. Gentle manipulation of the follicle may assist separation of oocytes from the follicle wall.
- 8 If flushing is required, the bung should be removed from the collection tube. A syringe containing flushing medium can be connected to the Luer connector, and flushing of the follicle performed. Upon completion of the flushing, the syringe should be removed and the bung reattached to the collection tube. Aspiration of the follicle should then be repeated.
- 9 Repeat for all follicles within the ovary. To minimise trauma and possible bleeding keep the needle tip within the ovary between aspirations.
- 10 Remove the collection tube from the silicone bung when full. Replace with a new sterile collection tube.
- 11 When all the follicles in one ovary have been aspirated, repeat the procedure with the other ovary. When all accessible follicles have been emptied, aspirate any blood or fluid that is present within the pouch of Douglas.
- 12 Withdraw the transducer and needle, check the vagina for signs of bleeding and allow the patient to rest. If vaginal bleeding is experienced post collection, apply steady pressure with a sponge and forceps or use an absorbable suture if bleeding persists.
- 13 Dispose of the needle set in accordance with local medical hazardous waste practices.

The above are only suggested instructions for use. The final decision on the technique used is the responsibility of the clinician in charge.

Wallace, Wallace design mark and Smiths Medical design mark are trademarks of Smiths Medical. The symbol © indicates the trademark is registered in the U.S. Patent and Trademark Office and certain other countries. All other names and marks mentioned are the trade names, trademarks or service marks of their respective owners. © 2015 Smiths Medical. All rights reserved.

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en

Single use.

Do Not Reuse: Medical devices require specific material characteristics to perform as intended. These characteristics have been verified for single use only. Any attempt to re-process the device for subsequent re-use may adversely affect the integrity of the device or lead to deterioration in performance. Sterile unless unit container is opened or damaged.

Destroy after single use.
Do not resterilize.

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Wallace, designmerket Wallace og designmerket Smiths Medical er varemerker som tilhører Smiths Medical. Symbolet ® angir at varemerket er registrert ved U.S. Patent and Trademark Office og i visse andre land. Alle andre navn og merker som er nevnt, er varemerker, handelsnavn og servicemerker som tilhører de respektive eiere.
Wallace, znak towarowy Wallace oraz znak towarowy Smiths Medical są znakami handlowymi Smiths Medical. Symbol ® oznacza, że znak towarowy jest zarejestrowany w Urzędzie Patentowym Stanów Zjednoczonych (U.S. Patent and Trademark Office) i w urzędach patentowych innych krajów. Wszystkie inne wymienione nazwy i znaki są nazwami towarowymi, znakami towarowymi lub nazwami usług ich właścicieli.