Dual Lumen Oocyte Recovery Set

en

INSTRUCTIONS FOR USE

fr

MODE D'EMPLOI

de

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it

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es

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pt

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nl

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sv

BRUKSANVISNING

da

BRUGSANVISNING

no

BRUKSANVISNING

fi

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el

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Tel: +44 (0)1303 260551

Australian Representative:

Smiths Medical Australasia Pty. Ltd.

Brisbane, QLD 4113, Australia.

Tel: +61 (0)7 3340 1300

www.smiths-medical.com

smiths medical

Wallace®

Dual Lumen Oocyte Recovery Set

INSTRUCTIONS FOR USE

These instructions contain important information for safe use of the product. Read the entire contents of these Instructions For Use, including Warnings and Cautions, before using this product. Failure to properly follow warnings, cautions and instructions could result in death or serious injury to the patient and/or clinician.

NOTE: DISTRIBUTE THIS INSERT SHEET TO ALL PRODUCT LOCATIONS.

Wallace® Dual Lumen Oocyte Recovery Sets are single-use devices provided sterile for ultrasound-guided, transvaginal collection of oocytes from the ovarian follicles.

WARNING:

- **READ THE COMPLETE DIRECTIONS FOR USE BEFORE USE. FAILURE TO PROPERLY FOLLOW THE INSTRUCTIONS, WARNINGS AND CAUTIONS MAY LEAD TO SERIOUS SURGICAL CONSEQUENCES OR INJURY TO THE PATIENT.**
- **THESE PROCEDURES SHOULD ONLY BE PERFORMED BY PERSONS HAVING ADEQUATE TRAINING AND FAMILIARITY WITH THESE TECHNIQUES. CONSULT MEDICAL LITERATURE REGARDING TECHNIQUES, COMPLICATIONS AND HAZARDS PRIOR TO PERFORMANCE OF THESE PROCEDURES.**
- **TO BE USED BY, OR UNDER THE DIRECTION OF, QUALIFIED PERSONS IN LINE WITH LOCAL GUIDELINES GOVERNING *IN VITRO* FERTILISATION, IF APPLICABLE.**
- **THE SHARP NEEDLE MUST BE HANDLED CAREFULLY TO AVOID NEEDLE-STICK INJURY, WHICH MAY RESULT IN TISSUE DAMAGE OR CROSS-INFECTION.**
- **DISCARD IF PRODUCT OR PACKAGING IS DAMAGED.**

Description

These instructions are intended for use with sterile Wallace® Dual Lumen Oocyte Recovery Sets, product codes:

Code	Description
DNS1633-950	16g 33cm dual lumen needle with 95cm transfer tubing
DNS1633-750	16g 33cm dual lumen needle with 75cm transfer tubing
DNS1633-500	16g 33cm dual lumen needle with 50cm transfer tubing
DNS1733-950	17g 33cm dual lumen needle with 95cm transfer tubing
DNS1733-750	17g 33cm dual lumen needle with 75cm transfer tubing
DNS1733-500	17g 33cm dual lumen needle with 50cm transfer tubing

The devices consist of a dual lumen stainless steel needle with 1cm of echo marking at the distal end. Attached to the proximal end are a plastic hub and a tubing set. The tubing set consists of a transfer tube leading through a silicone bung, and a vacuum tube also passing through the bung and terminating in a Luer compatible connector. Attached to the plastic hub is a flushing tube, terminating in a clear Luer compatible connector.

A pump adaptor is supplied in the pack to connect the vacuum tube to a vacuum source (if required), and is compatible with elastic tubing internal diameter 4mm-10.5mm.

The devices are colour coded. Blue for 16g, Red for 17g.

Indications for use

Ultrasound-guided transvaginal collection of oocytes.

Contraindications

The Dual Lumen Oocyte Recovery Set should not be used:

- In the presence of or after recent pelvic inflammatory disease.
- For intra-fallopian procedures.

Instructions for use

- 1 Oocyte recovery should be performed under either general anaesthetic or analgesia.
- 2 Place the patient in the lithotomy position, prepare the vulva and vagina for full aseptic technique and cover the area in sterile drapes.
- 3 Cover the ultrasound transducer with a sterile sleeve with a small amount of

coupling gel covering the end of the probe within the sleeve. Introduce into the vagina and scan the pelvis to establish the position of the ovary, the development of the endometrium and the number and position of follicles to be aspirated.

- 4 Remove the oocyte recovery set from its packaging and connect a syringe or vacuum pump to the Luer compatible connector on the vacuum tube (use the pump adaptor if required) and a sterile collection tube to the silicone bung.
- 5 Connect a sterile syringe or other suitable delivery system containing flushing media to the Luer compatible connector at the end of the flushing tube.
- 6 Flush the transfer and flushing tubes/lumen with media. After flushing purge the transfer system of fluid and empty or discard the collection tube. Prime the flushing tube and flushing lumen with media.

WARNING

Test the system to ensure all connections are sound and flow can be achieved. Failure to test the system may result in failure to obtain the required sample and/or tissue damage.

- 7 Insert the needle into the needle guide of the ultrasound transducer, ensuring that there is smooth movement along its length.
- 8 Insert the ultrasound transducer and needle into the vagina and direct the probe towards the posterolateral aspect of the vaginal fornix. Manipulate the probe until the follicles appear in close proximity to the needle.

WARNING

Failure to use ultrasound to locate the follicles or incorrect needle guidance may result in failure to obtain the required sample and/or tissue damage.

- 9 Advance the needle into the nearest follicle. If the ovaries are mobile apply pressure suprapubically and enter the follicles. When the tip of the needle can be seen in the follicle, aspirate the contents using either the syringe or vacuum pump attached to the set. Gentle manipulation of the follicle may assist separation of oocytes from the follicle wall. An initial vacuum pump setting in the range 100-200mmHg is suggested, the vacuum source should be controlled to give the desired flow rate for each patient.
- 10 If flushing is required, inject media into the follicle using either the syringe or delivery system connected to the flushing tube, until the follicle returns to its original size. Aspiration of the follicle should then be repeated.
- 11 Repeat for all follicles within the ovary. To minimise trauma and possible bleeding keep the needle tip within the ovary between aspirations.
- 12 Remove the collection tube from the silicone bung when full (do not allow the collection tube to become so full that oocytes are aspirated into the vacuum source). Replace with a new sterile collection tube.
- 13 When all the follicles in one ovary have been aspirated, repeat the procedure for the other ovary. If necessary, once all accessible follicles have been emptied, aspirate any blood or fluid that is present within the pouch of Douglas.
- 14 Withdraw the transducer and needle, check the vagina for signs of bleeding and allow the patient to rest. If vaginal bleeding is experienced post collection, apply steady pressure with a sponge and forceps or use an absorbable suture if bleeding persists.
- 15 Dispose of the needle set in accordance with local medical hazardous waste practices.

The above are only suggested instructions for use. The final decision on the technique used is the responsibility of the clinician in charge.

Wallace, Wallace design mark and Smiths Medical design mark are trademarks of Smiths Medical. The symbol ® indicates the trademark is registered in the U.S. Patent and Trademark Office and certain other countries. All other names and marks mentioned are the trade names, trademarks or service marks of their respective owners. © 2012 Smiths Medical. All rights reserved.

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en**Single use.**

Do Not Reuse: Medical devices require specific material characteristics to perform as intended. These characteristics have been verified for single use only. Any attempt to re-process the device for subsequent re-use may adversely affect the integrity of the device or lead to deterioration in performance.

Sterile unless unit container is opened or damaged. Destroy after single use. Do not resterilize.

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Wallace, Designmerket Wallace og designmerket Smiths Medical er varemerker som tilhører Smiths Medical. Symbolet ® angir at varemerket er registrert ved U.S. Patent and Trademark Office og i visse andre land. Alle andre navn og merker som er nevnt, er varemerker, handelsnavn og servicemerker som tilhører de respektive eiere.
Wallace, Wallace a Smiths Medical sú ochranné známky Smiths Medical. Symbol® označuje, že ochranná známka je registrovaná na úrade U.S. Patent and Trademark Office a v istých ďalších štátoch. Všetky ďalšie uvádzané názvy a známky sú obchodnými názvami, ochrannými značkami alebo servisnými značkami svojich príslušných vlastníkov.
Wallace, Znak towarowy Wallace oraz znak towarowy Smiths Medical są znakami handlowymi Smiths Medical. Symbol® oznacza, że znak towarowy jest zarejestrowany w Urzędzie Patentowym Stanów Zjednoczonych (U.S. Patent and Trademark Office) i w urzędach patentowych innych krajów. Wszystkie inne wymienione nazwy i znaki są nazwami towarowymi, znakami towarowymi lub nazwami usług ich właścicieli.