Intra-uterine Insemination Catheters

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www.smiths-medical.com

CCN. 004/006/008A 07/14

Wallace® Intra-uterine Insemination Catheters

INSTRUCTIONS FOR USE

Wallace® intra-uterine insemination catheters are single-use devices provided sterile for the introduction of washed spermatozoa into the uterine cavity.

These instructions contain important information for safe use of the product. Read the entire contents of these Instructions For Use, including Warnings and Cautions, before using this product. Failure to properly follow warnings, cautions and instructions could result in death or serious injury to the patient and/or clinician.

NOTE: DISTRIBUTE THIS INSERT SHEET TO ALL PRODUCT LOCATIONS.

WARNING:

- **THESE PROCEDURES SHOULD ONLY BE PERFORMED BY PERSONS HAVING ADEQUATE TRAINING AND FAMILIARITY WITH THESE TECHNIQUES. CONSULT MEDICAL LITERATURE REGARDING TECHNIQUES, COMPLICATIONS AND HAZARDS PRIOR TO PERFORMANCE OF THESE PROCEDURES.**
- **TO BE USED BY, OR UNDER THE DIRECTION OF, QUALIFIED PERSONS IN LINE WITH LOCAL GUIDELINES GOVERNING IN VITRO FERTILISATION, IF APPLICABLE.**
- **DISCARD IF PRODUCT OR PACKAGING IS DAMAGED.**

Description

These instructions are intended for use with sterile Wallace® intra-uterine insemination catheters, product code AIC18.

The devices consist of a flexible inner catheter and an outer sheath which is attached to the inner catheter by a Luer compatible connector.

The inner catheter is 16g with a lumen that is uniform throughout its length, and capable of containing 0.2ml of fluid. The catheter has 2 side eyes located at the smooth, closed distal tip.

The inner catheter protrudes from the outer sheath by 5cm, and the outer sheath has a series of 1cm graduations at the distal end.

Indications for use

Introduction of washed spermatozoa into the uterine cavity. (Intra-uterine artificial insemination)

Contraindications

The catheter should not be used:

- In the presence of chronic cervical infection
- In the presence of or after recent pelvic inflammatory disease
- For intra-fallopian procedures

Instructions for Use

- 1 Place patient in the lithotomy, dorsal, or left lateral position. The cervix should be exposed with a speculum and gently swabbed with cotton wool moistened with normal saline or medium.

When using transabdominal ultrasound follow normal procedures such that an optimal image of the uterine cavity is achieved.

- 2 With the plunger fully depressed firmly attach a syringe (preferably 3ml or 5ml) to the inner catheter.
- 3 Draw 1ml of air into the syringe, then place the distal tip of the catheter into the specimen of washed spermatozoa so that both lateral ports are immersed. Aspirate the desired quantity of specimen into the catheter.

WARNING

Always use washed spermatozoa when performing intra-uterine insemination. The introduction of unwashed spermatozoa into the uterus will almost always result in severe adverse reaction which may include anaphylactic shock. Refer to medical literature for methods of preparing spermatozoa for intra-uterine insemination before performing this procedure.^[1]

- 4 With the hubs locked in position, advance the catheter so that the inner catheter passes through the external and internal os, into the mid-uterine cavity. It may be necessary to twist the catheter as it negotiates the cervical canal. When using ultrasound manipulate the probe to obtain the best image of the inner catheter.

WARNING

The catheter should never be forced against digitally felt resistance while inside the uterine cavity, as forcing the catheter may result in damage to the endometrial tissue and bleeding.

The outer sheath should not be advanced further than the internal os, and should certainly never enter the uterine cavity, as this may result in damage to the endometrial tissue and bleeding.

The inner catheter tip should not be pulled back further than the end of the outer sheath, as this may result in damage to the endometrial tissue and bleeding.

- 5 If slight resistance is experienced, withdraw the catheter and advance the outer sheath until only the tip of the inner catheter is visible. Using thumb and forefinger curve the outer sheath to complement the patient's anatomy and pass the catheter in this position to the internal os. In the rare event where the cervical canal cannot be negotiated by this method, stylets are available to assist insertion. Use code no 18165T. Advance the inner catheter into the mid-uterine cavity.
- 6 Place tip of catheter approximately 1cm from the fundus and slowly inject the washed spermatozoa and 1ml of air in the syringe to ensure full expulsion of the specimen. Should there be any doubt about the complete expulsion, remove the syringe from the catheter; reload with air; reattach the syringe and slowly inject the air to expel any remaining specimen.
- 7 Leave time for the uterine environment to equilibrate and slowly remove the catheter.
- 8 Remove the speculum and allow the patient to rest.
- 9 Dispose of the catheter in accordance with local medical hazardous waste practices.

The above are only suggested instructions for use. The final decision on the technique used is the responsibility of the clinician in charge.

References

- [1] Lopata A, Patullo MJ, Chang A, James B 1976
A method for collecting motile spermatozoa from human semen. *Fertil Steril* 27:677

Wallace, Wallace design mark and Smiths Medical design mark are trademarks of Smiths Medical. The symbol ® indicates the trademark is registered in the U.S. Patent and Trademark Office and certain other countries. All other names and marks mentioned are the trade names, trademarks or service marks of their respective owners.

© 2014 Smiths Medical. All rights reserved.

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Wallace, designmerket Wallace og designmerket Smiths Medical er varemerker som tilhører Smiths Medical. Symbolet © angir at varemerket er registrert ved U.S. Patent and Trademark Office og i visse andre land. Alle andre navn og merker som er nevnt, er varemerker, handelsnavn og servicemerker som tilhører de respektive eiere.
Wallace, znak towarowy Wallace oraz znak towarowy Smiths Medical są znakami handlowymi Smiths Medical. Symbol® oznacza, że znak towarowy jest zarejestrowany w Urzędzie Patentowym Stanów Zjednoczonych (U.S. Patent and Trademark Office) i w urzędach patentowych innych krajów. Wszystkie inne wymienione nazwy i znaki są nazwami towarowymi, znakami towarowymi lub nazwami usług ich właścicieli.