



ORIGIO® VITRIFICATION

# VitriFit™ - Open Carrier

A cool design for optimal storage



**VitriFit**

Now available in  
packs of 5 and  
8 colors

CE  
0543


**CooperSurgical**  
Fertility and Genomic Solutions

## Maximize efficiency using minimal storage space with VitriFit


VitriFit is an open system vitrification carrier for storage of oocytes, embryos and blastocysts. It works for all stages and with all media kits.

The hexagonal shape allows VitriFit carriers to be packed tighter during storage than other popular vitrification carriers on the market, thereby improving storage efficiency by up to 50%.

The handle and cap are made from the same material and are designed for a perfect fit at all temperatures. The loading area is curved for specimen protection and to ease loading and removal of excess medium to reach minimal volume.



## Survival rates for specimen obtained with VitriFit compared to popular open carrier alternatives



### Reference Numbers



Individually Wrapped (20 pouches per box)	Pack of 5 (10 pouches per box)
42782001 - Clear	42785005 - Clear
42792001 - Blue	42795005 - Blue
42802001 - Green	42805005 - Green
42812001 - Yellow	42815005 - Yellow
42822001 - Lime	42825005 - Lime
42842001 - Purple	42845005 - Purple
42852001 - Orange	42855005 - Orange
42862001 - Pink	42865005 - Pink