


Reassuring Quality

TPC Micropipettes - Designed by Embryologists for Embryologists



When you need the tool that perfectly fits the task at hand, choose TPC.

CooperSurgical®

Fertility and Genomic Solutions

ICSI Micropipettes

Type	Reference no.	Toe length (tip to bend distance)	Bend angle	ID	OD	Bevel angle
Injection (no spike)	LICR-BA21TA35	0.6 mm	35	4.5-5.5	6.0-7.0	21°
Injection	LICR-ID4	1 mm	30	4.0-4.5	5.0-6.0	24°
Injection	LICR-ID4TA30	0.6 mm	30	4.0-4.5	5.0-6.0	24°
Injection (no spike)	LICR-OD5	1 mm	30	4.0-4.5	5.0-6.0	20°
Injection	LICR-TA20	1 mm	20	4.5-5.5	6.0-7.0	30°
Injection	LICR-TA35	1 mm	35	4.5-5.5	6.0-7.0	30°
Injection	LICR-TL0.5	0.5 mm	30	4.5-5.5	6.0-7.0	30°
Injection	LISR	1 mm	30	4.5-5.5	6.0-7.0	30°

Holding Micropipettes

Type	Reference no.	Toe length (tip to bend distance)	Bend angle	ID	OD
Holding	LHC-ID15	1 mm	30	14-18	110-130
Holding	LHC-ID30	1 mm	30	28-35	110-130
Holding	LHC-ID30TA35	1 mm	35	28-35	110-130
Holding	LHC-OD100ID20	1 mm	30	18-23	90-110
Holding	LHC-OD100TA20	1 mm	20	18-23	90-110
Holding	LHC-TA20	1 mm	20	20-27	110-130
Holding	LHC-TA35	1 mm	35	20-27	110-130
Holding	LHS	1 mm	30	20-27	110-130

Biopsy Micropipettes

Type	Reference no.	Toe length (tip to bend distance)	Bend angle	ID	OD
Biopsy	LBC-OD27BA90	1 mm	30	19-23	25-29
Biopsy	LBC-OD30BA90	1 mm	30	21-25	28-32
Biopsy	LBC-OD32TA35	1 mm	35	23-27	30-34
Biopsy	LBC-OD35BA90	1 mm	30	25-29	33-37



*Day-to-day delivery may vary according to geographical location